



Cirrus 1700 & 1500 Installation Guide

The Cirrus bath is supplied with the following:

Front panel, twin waste kit, wall fixing brackets

Optional items: End panel, bath taps, antiscald valve

IMPORTANT NOTES

The bath is fitted with 6 supporting feet which are to be adjusted to accommodate the front panel. The feet should be adjusted to fully support the bath floor and avoid any distortion in use.

Twin wastes are supplied to increase drainage and to provide security if one waste should fail.

If fitting in an alcove, ensure the bath goes freely into position. Sideways pressure can distort the doorway.

The bath is manufactured in glass fibre with a gel-coat surface. The bath edges are more brittle than acrylic baths.

STANDARD FITTING GUIDE

- 1) Fit the wall fixing brackets to bath timber frames at the rear of the bath.
- 2) Drill the bath and fit the taps.
- 3) Fit the wastes, overflow and waste traps.
- 4) Place the bath into position & adjust to height (approx. 550mm). Use the front panel as a guide. It is advisable to set the feet onto battens on to the floor to avoid over extending the legs.
- 5) Level the bath in all planes and lock the feet.
- 6) Mark wall fixings and plumbing runs.
- 7) Remove the bath, drill wall fixings and run plumbing into position.
- 8) Re-position bath and affix to wall and to floor. *Ensure the 2 bath feet under the doorway are adjusted to support the foot well as well as the doorway integrity.*
- 9) Connect plumbing. The wastes may be run independently or joined by use of a tee. Refer to fitting guides for the thermostatic valve.
- 10) Test plumbing for leaks.
- 11) Fill the bath to maximum and leave for about 30 minutes to test the door seal.

- 12) Whilst the bath seal is tested, position and fit the front panel (and end panel if required).
The front panel fits into clips under the top rim. The panel base will require a batten. The end panel will require battens along floor and on the wall.
- 13) Secure the front and end panels using screws and caps. Do not use silicone sealant as this prevents future service access.
- 14) Finally drain bath and **leave the door open** to prevent compression of the seal.
- 15) Instruct the customer on the use of the bath.

◆ Refer to User Guide and leave copy with user.

◆ **Do not leave door in locked position on completion. The door should be left unlocked when not in use in order to preserve the seal.**

September 2020